



Simply connect and get started!

YOUNG PROXEO TWIST

wireless handpiece with foot pedal

The Young Proxeo Twist Cordless enables rotary polishing without any restrictions – in more ways than one. **Complete cordless freedom of movement** ensures maximum comfort during use. The handpiece has a cordless foot control and is compatible with two different cup systems: Young polishing cups with screw connection or hygienic disposable contra-angles.

**Recommended
rotation speeds**
1.250 - 2.500 rpm



Modular
Choose from various
cup systems



Versatile
easy to move from
room to room

Well balanced
Reduces hand fatigue

 **Cordless and
lightweight**

Reliable
Li-ion battery
lasts all day



The Benefits:

Ideal rotation speeds protect the tooth from overheating and prevents abrasion. It also ensures a comfortable treatment and reduces polishing paste splatter.

The **wireless foot control** allows maximum freedom of movement and quick position changes around the patient. Easy to operate thanks to WiFi connection. The speed can be adjusted continuously from 0 to 2,500 rpm.


Ergonomic treatments are made possible by **freedom of movement** working. With cordless freedom of movement the treatment position can be adjusted at any time and ensures maximum comfort during use. The cordless handpiece also enables a back-friendly working posture. **Also suitable for left-handed users.**

Item number	Description
YI-PL-64SET	Young Proxeo Twist Cordless Set, consisting of drive handpiece, handpiece holder, charger with adapter, foot pedal and screw-type attachment YI-PL-64
YI-PL-44SET	Young Proxeo Twist Cordless Set, consisting of drive handpiece, handpiece holder, charger with adapter, foot pedal and attachment for disposable contra-angles YI-PL-44
YI-PL-64	Young Proxeo Twist Cordless Contra-angle Handpiece Attachment for Screw-Type Polishing Cups
YI-PL-44	Young Proxeo Twist Cordless Handpiece Attachment for Disposable Contra-Angles, Doriot System

The following products are included in the kits (YI-PL-64SET & YI-PL-44SET)



Discover all the benefits now



YI-D-PTCH/1125